bio-techne[®]

CHARACTERIZATION OF MAURICE[™] CE-SDS PLUS FOR USP <129> SUITABILITY

Pouya Amrollahi, Will McElroy, and Chris Heger

ProteinSimple, a Bio-Techne Brand, 3001 Orchard Parkway, San Jose, CA 95134, USA

INTRODUCTION

Since the first therapeutic monoclonal antibody (mAb) was commercialized in the mid-80's, close to 100 therapeutic mAb products (accounting for around a quarter of all biotech drugs) have hit the market; making it a \$125 billion industry that targets critical pathological health conditions – including but not limited to products for antitumor, antiviral, and antiplatelet therapies. From early-stage process development to batch lot release testing, the efficacy, safety, identity, stability, and purity of therapeutic mAb products throughout their shelf life are of crucial importance. Capillary electrophoresis sodium dodecyl sulfate (CE-SDS) has become the gold standard technique for the quality-control of therapeutic mAbs and proteins due to its ease of implementation, robustness, and reproducibility, replacing the more traditional and labor-intensive technique such as SDS-PAGE gel. Successful CE-SDS method development, under both reducing and nonreducing conditions, aims to reduce assay-associated impurities, fragmentations, and aggregations.

Here, we have used the monoclonal IgG System Suitability Reference Standard developed by U.S. Pharmacopeia (USP) to assess the rigor and robustness of an optimized Maurice[™] CE-SDS PLUS method compared to the recommended USP protocol provided in monograph <129>. The optimization leveraged Design of Experiments (DOE) to optimize key components in sample preparation, denaturing conditions, and sample injection. The results show that the optimized methods: (1) cause less fragmentation compared to the USP <129> method, (2) are not susceptible to sample injection variations that might differ between instruments, and (3) provide comparable data to the USP <129> monograph for mAbs.

MATERIALS & METHODS

	aunto	REAGENT/KIT	VENDOR	PART #
proteinsimple		Maurice [™] CE-SDS <i>PLUS</i> Application Kit	ProteinSimple	PS-MAK03-S
Maurice		Monoclonal IgG System Suitability	U.S. Pharmacopeia	1445550
-	CE-SDS PLUS Cartridge	β-mercaptoethanol (β-ME)	Sigma-Aldrich	M-3148
Maurice™		Iodoacetamide (IAM)	Sigma-Aldrich	16125

MAURICE[™] CE-SDS METHOD

Lyophilized monoclonal IgG system suitability was reconstituted with CE-SDS PLUS buffer to a concentration of 1mg/mL in a final volume of 50µL, mixed with 2µL of reconstituted 25X Internal Standard (included in the kit), and either 2.5µL of 14.2M β-ME (reduced IgG samples) or 2.5µL of 20mM IAM (non-reduced IgG samples) before heat denaturation (10min at 70°C for reduced IgG samples and 5min at 65°C for nonreduced IgG samples) in a thermocycler. Samples were kept on ice for 5min before transferring to a 96-well plate, then centrifuged for 10min at 1000×g, and inserted in Maurice[™]. Batch reagents (included in the kit) were loaded into Maurice[™] based on the application guide. Samples were injected for 20sec at 4600V and separated for 25min (reduced IgG) or 35min (non-reduced IgG) at 5750V using the CE-SDS PLUS cartridge.

DESIGN OF EXPERIMENT AND DATA ANALYSIS

Suitable statistical methods, Box-Behnken and D-optimal designs, were applied using the JMP® software for design of experiment. Sample concentration, sodium dodecyl sulfate (SDS) concentration, β -ME and IAM concentration, denaturation condition for reduced and non-reduced methods, separation time, and voltage were optimized to ensure complete reaction, while minimizing method-induced fragmentation. Assay robustness was tested through 48-injection batches and the methods were compared to the recommended U.S. Pharmacopeia protocol <129>. All data were analyzed with Compass for iCE software and JMP[®].

RESULTS

We started the method optimization by first comparing commonly used CE-SDS sample buffers. We used SCIEX[™] sample buffer, CE-SDS buffer, and CE-SDS PLUS To study the robustness and reproducibility of our optimized assays, we ran 48injection batch experiments. We evaluated robustness by looking at the relative buffer for this purpose. CE-SDS PLUS buffer outperformed the other two options in both reduced and non-reduced experiments with higher injection efficiency standard deviation (%RSD) of peak area of heavy chain (%RSD: %2.5) and intact (Figure 1). We further evaluated the robustness of methods by probing sample (Sample concentration: 0.6-1.4mg/mL and Buffer SDS concentration: 0.5-1.5X of CE-SDS PLUS buffer) and instrument (Separation time: 10-30 min and Injection voltage: 4500-5500V) variables. While sample concentration was linearly correlated with peaks (%RSD: %2.7) in reduced and non-reduced methods, respectively (Figure 5). fragmentation as expected, we saw strong assay performance through other parameters; showing method suitability across instruments and SDS concentrations (Figure 2).



Figure 1. Evaluating different sample buffers. Compared to the SCIEX[™] and CE-SDS buffer, CE-SDS PLUS buffer helps injecting the species of interest better in each method; resulting in the highest area for intact and heavy chain peak.

Reduced Method

A successful reduced CE-SDS method is achieved when denaturation is A non-reduced CE-SDS method is used for evaluating sample purity, measuring Finally, we compared the optimized CE-SDS *PLUS* methods with the analytical completed with minimal induced fragmentation and background peaks. Higher intact product, fragments, covalently bound aggregates, and non-product procedure for recombinant therapeutic mAb described in U.S. Pharmacopeia injection voltage and elevated temperatures improve resolution. We sought to related impurities. mAbs are incubated at high temperatures for denaturation <129> protocol. We ran three replicates for each method (total of 12 samples) in optimize the reducing agent (here β -ME) concentration, denaturation time, and (accelerates SDS binding) and alkylation; however, it can induce mAb the same run where each sample was injected four times. Results showed that denaturation temperature. We probed 0.3-1.1M of β -ME and denatured the fragmentation. We sought to optimize alkylating agent (here IAM) the optimized CE-SDS PLUS methods match the sensitivity and reproducibility of samples at 65-75°C for 5-15min. Samples treated at a temperature lower than concentration, as well as denaturation time and temperature, by minimizing the USP <129> standard protocol. There were no significant differences between the 70°C showed extra peaks associated with incompletely reduced IgG, that ran peak area of fragments. Our results showed that the most significant factor was amount of incomplete reduction in reduce methods (%CE-SDS PLUS: 0.98±0.09 later than heavy chain peak (Figure 3). Our results showed that increasing the IAM concentration. Moreover, denaturing samples at lower temperatures and and %USP <129>: 1.00±0.04). However, having lower method-induced amount of reducing agent and denaturation time did not change the profile; for shorter time led to a \sim 15% reduction in the peak area of fragments (Figure fragmentation, the non-reduced CE-SDS *PLUS* method showed significant therefore, we selected the following optimized CE-SDS PLUS method conditions: 4); therefore, we selected the following optimized CE-SDS PLUS method conditions: improvements (*P*-value < 0.0001) over the USP non-reduced method (Figure 6).



biotechne brands

Maurice[™] CE-SDS *PLUS* Method Optimization

Figure 2. Evaluating sample and instrument conditions. Sample concentration is linearly correlated with fragmentation However, our method is not susceptible to changes of %SDS in CE-SDS PLUS buffer, as well as instrument parameters.

Non-reduced Method







Assay robustness Performing 48 injection batch experiments, using both non-reduced and reduced methods showed strong assay robustness with low degree of njection-to-injection deviation during long hours of experiment (48 hours). The slight drop after every 12 injection is due to the cartridge clean-up step in our method.

Comparison to USP <129> Protocol



Figure 6. Comparing CE-SDS PLUS method to USP <129> standard protocol. Statistical analysis (Mann-Whitney test) shows that non-reduced CE-SDS PLUS method induces significantly less fragmentation compared to USP <129> method; while both reduced methods perform equally acceptable in minimizing incomplete reduction (ns: not significant).

CONCLUSIONS

- Maurice[™] CE-SDS *PLUS* methods developed here show robust performance over 48-hour-long CE-SDS experiments.
- The optimized reduced method matches the sensitivity of USP <129> protocol, while the non-reduced method significantly outperforms it.



CONTACT: pouya.amrollahi@bio-techne.com

bio-techne.com © 2020 Bio-Techne[®]. All rights reserved.