



USAID
FROM THE AMERICAN PEOPLE

Promoting the
QUALITY OF MEDICINES Plus

Chlorhexidine Digluconate (7.1%) Gel Job Aid to Assist with Dossier Preparation

January 2021



Contact Information

Promoting the Quality of Medicines Plus Program
United States Pharmacopeia
12601 Twinbrook Parkway
Rockville, MD 20852 USA
Tel: +1-301-816-8166
Fax: +1-301-816-8374
Email: PQMplus@USP.org

This document is made possible by the generous support of the American people through the U.S. Agency for International Development (USAID) Cooperative Agreement No. AID-7200AA19CA00025. The contents are the responsibility of U.S. Pharmacopeial Convention (USP) and do not necessarily About PQM+

The Promoting the Quality of Medicines Plus (PQM+) Program is a five-year cooperative agreement between USAID and USP to sustainably strengthen medical product quality assurance systems in low- and middle-income countries. The program works to improve medical product quality through cross-sectoral and systems strengthening approaches and the application of international quality assurance standards across the pharmaceutical system. By sharing scientific expertise and providing technical support and leadership, PQM+ helps create resilient and robust local health systems that address diseases such as HIV/AIDS, tuberculosis, malaria, and neglected tropical diseases, as well as improve maternal, newborn, and child health.

Suggested Citation

This document may be reproduced if credit is given to PQM+. Please use the following citation:

PQM+. 2020. Chlorhexidine digluconate (7.1%) Gel Dossier Aid. Submitted to the U.S. Agency for International Development by the PQM+ Program. Rockville, MD: U.S. Pharmacopeial Convention

Acknowledgments

The original report (Chlorhexidine digluconate (7.1%) gel Technology Transfer Report – May 2019), from which information in this job aid was compiled, was written by USP in collaboration with GSK. USP gratefully acknowledges the preparation of this job aid by Frederick Meadows, Senior Technical Advisor, Product Supply Management and CMC, PQM+ under the direction of Lawrence Evans, Senior Director, Technical PQM+.

USP is also grateful to the U.S. Agency for International Development (USAID), which funded this work through the Promoting the Quality of Medicines Plus (PQM+) program. Specifically, USP is thankful for the guidance and oversight provided by:

- Alison Collins, Health Systems Advisor
- Lisa Ludeman, Senior Pharmaceutical Management Advisor
- Tobey Busch, Senior Pharmaceutical Management Advisor
- Helen Petach, Senior MCH Science Advisor, USAID.

Contents

Acronyms.....	Er
ror! Bookmark not defined.	
Aid for CTD Module 2.3.A.1 – Facilities and Equipment.....	5
Aid for CTD Module 3.2.S.1.3 – API Physicochemical Attributes.....	5
Aid for CTD Module 3.2.S.3.2 – Impurities.....	5
Aid for CTD Module 3.2.P.2.3 – Manufacturing Process Development.....	8
Aid for CTD Module 3.2.P.3.2 – API and CHX Gel Formulation.....	8
Aid for CTD Module 3.2.P.3.3 – Impurity Control Strategy.....	8
Aid for CTD Module 3.2.P.3.4 – Analytical and Specifications for CHX Gel.....	9
Aid for CTD Module 3.2.P.7 – Manufacturing and Packaging.....	10
References.....	11

Acronyms

4CA	4-Chloroaniline
API	active pharmaceutical ingredients
CHX	Chlorhexidine
CTD	common technical document
CQA	critical quality attribute
HPLC	high performance liquid chromatography
Ph. Eur.	European Pharmacopeia
PLC	programmable logic controller
PQM+	Promoting the Quality of Medicines Plus
USAID	U.S. Agency for International Development
USP	U.S. Pharmacopeial Convention

This document is a job aid to enable quick access to technical information by those preparing dossiers for Chlorhexidine (CHX) digluconate (7.1%) gel. The information consists of facilities, equipment, safety recommendations, and analytical data for production. The information herein is adapted from the Promoting Quality Medicines Program's Chlorhexidine digluconate (7.1%) gel Product Information Report (May 2019).

Aid for CTD Module 2.3.A.1 – Facilities and Equipment

Table 1. Facilities for Processing CHX Gel

Parameter	Comments
Permitted Daily Exposure Equipment and Facility Controls	<ul style="list-style-type: none"> • Consult toxicological information and regulatory authorities.
Compounding CHX Gel	<ul style="list-style-type: none"> • Vessel – A) suitably sized and jacketed for heating and cooling; B) sealable stainless steel; C) PLC control system, e.g., Siemens PCS7 controller • Variable speed, bidirectional homogenizer • Variable speed agitator • Vacuum and compressed air for material transfer • Stainless steel conical vessels for intermediate bulk storage (200–300 Kg capacity) for CHX gel
CHX Gel Bulk Product Equipment	<ul style="list-style-type: none"> • Proper geometry for 1) uniform mixing across vessel, e.g., spherical or conical that avoids dead spaces; and 2) good product recovery • Allows vacuum deaeration with controlled addition of liquids and solids • Mixing system for uniform mixing without vortices • External heating and cooling jacket

Aid for CTD Module 3.2.S.1.3 – API Physicochemical Attributes

Table 2. API Physicochemical Properties

Test	Comments
Appearance	<ul style="list-style-type: none"> • CHX digluconate solution, 20% w/v is a clear to pale yellowish liquid conforming to the European Pharmacopeia (Ph. Eur.) Monograph for description.
Identity and Assay	<ul style="list-style-type: none"> • Identity and Assay by ultraviolet and high performance liquid chromatography (HPLC) • Identity and Assay by HPLC per Ph. Eur. Monograph for description

pH	<ul style="list-style-type: none"> CHX digluconate active pharmaceutical ingredients (API) solution should be controlled within the range of 5.5–7.0 to minimize impurity formation, e.g., 4CA and others.
4CA (most critical impurity)	<ul style="list-style-type: none"> Impurity supplier specification \leq 20 ppm in CHX digluconate solution (20% w/v) Stays below 150 ppm \leq 6 months at Zone II (25°C/60% RH) Should be controlled within the range of 5.5–7.0 Impacts 4CA and other impurity formation and should be controlled for stability 4CA formed by pH and thermal stress 4CA is a CQA – genotoxic and carcinogenic
Other Impurities	<ul style="list-style-type: none"> Occur via hydrolysis, condensation, and other mechanisms Controlled in accordance with the levels permitted in the Ph. Eur. Monograph
Thermal Stability	<ul style="list-style-type: none"> Stable for 6 months under Zone II climatic conditions (25°C/60% RH)

Aid for CTD Module 3.2.S.3.2 – Impurities

Table 3. CHX Impurities

CQA Parameter	Release Specification (acceptance criteria)	Shelf-life Specifications (acceptance criteria)
4CA (Impurity P)	NMT 0.08%	NMT 0.35%
1-(6-aminohexyl)-5-(4-chlorophenyl) biguanide (Impurity G)	NMT 0.3%	NMT 0.6%
1-[6-(carbamimidoylamino)hexyl]-5-(4-chlorophenyl) biguanide (Impurity N)	NMT 1.0%	NMT 2.0%
1-(4-chlorophenyl)-5[6-[[4-[(4-chlorophenyl) amino]-6[(1S,2R,3R,4R)-1,2,3,4,5-pentahydroxypentyl]-1,3,5-triazin-2-yl] amino] hexyl] biguanide (Impurity J)	NMT 0.4%	NMT 0.7%
1,1-[iminobis(carbonimidoylimino)hexane-6,1-diyl] bis [5-(4-chlorophenyl) biguanide] (Impurity H)	NMT 0.5%	NMT 0.5%
1-(4-chlorophenyl)-5[6[cyanocarbamidoyl]amino] hexyl] biguanide (Impurity A)	NMT 0.4%	NMT 0.4%
Sum of Impurity (I+O), 5-(2-chlorophenyl)-5'-(4-chlorophenyl)-1&1'-(hexane-1,6-diyl) dibiguanide	NMT 0.4%	NMT 0.4%

N-(4-chlorophenyl)-N'[[6[[[(4-chlorophenyl) carbamimidoyl] amino] hexyl] carbamimidoyl] urea (Impurity K)	NMT 0.4%	NMT 0.7%
N-(4-chlorophenyl) urea (Impurity F)	NMT 0.2%	NMT 0.3%
Any other specified impurities*	NMT 0.2%	NMT 0.2%
Any other unspecified impurities	NMT 0.1%	NMT 0.2%
Total impurities	NMT 2.0%	NMT 4.0%
Microbial content	Complies with harmonized pharmacopeia	
Total aerobic microbial count	≤10 ² cfu/g	≤10 ² cfu/g
Total combined yeasts/mold count	≤10 ¹ cfu/g	≤10 ¹ cfu/g
Staphylococcus aureus	Absent in 1g	Absent in 1g
Pseudomonas aeruginosa	Absent in 1g	Absent in 1g

* Any other specified impurities include:

- 1-[6-[[[(4-chlorophenyl) carbamimidoyl] carbamimidoyl] amino] hexyl] carbamimidoyl] urea (Impurity B)
- (5R,6S)-2-[(4-chlorophenyl) amino]-5-hydroxy-6-[(1R,2R)-1,2,3-trihydroxypropyl]-5,6-dihydro-4H-1,3-oxazin-4-one (Impurity L)
- Impurity Q

Table 4. CHX Elemental Impurities

Substance	Specification
Elemental Impurities (per ICH Q3D)	<ul style="list-style-type: none"> • Risk assessment should be conducted to include: A) input API; B) excipients; C) processing aids (water); D) manufacturing equipment; and E) contain closure systems. The oral route guidance should be used for CHX gel. • Risk assessment focuses on Class 1 (cadmium, lead, arsenic, and mercury); Class 2A (cobalt, vanadium, and nickel) elemental impurities. • Likelihood of reaching final dose should be deemed to be low. • Very low levels of relevant elements, <10% of permitted daily exposure limit, can be approved in final product.

Aid for CTD Module 3.2.P.2.3 – Manufacturing Process Development

Table 5. CHX Manufacturing, Filling and Sealing

Unit Operation	Comments
Gel Manufacture	<ul style="list-style-type: none"> Disperse and hydrate guar gum with heating in an aqueous solution of sodium acetate trihydrate. Add in CHX digluconate API solution. Mix.
Filling/Sealing	<ul style="list-style-type: none"> Fills into aluminum foil laminate sachets by automated filling/sealing equipment. Target fill weight for sachets is 3.3g (range 3.0g–3.6g).

Aid for CTD Module 3.2.P.3.2 – API and CHX Gel Formulation

Table 6. CHX Gel Formulation Composition

Substance	Content
API (salt forms)	<ul style="list-style-type: none"> dihydrochloride (0.2% aqueous solubility) diacetate (2% aqueous solubility) digluconate (highest solubility)
CHX digluconate (7.1%) gel	<ul style="list-style-type: none"> 20% w/v aqueous solution (37.81% w/w) Consists of three excipients (guar gum (1.40 % w/w) as a thickener, sodium acetate trihydrate (0.10 % w/w) as a pH modifier, QS with purified water as a solvent) 4% w/w equivalent of CHX as the free base Single-use 3g dose in a foil laminate sachet

Aid for CTD Module 3.2.P.3.3 – Impurity Control Strategy

Table 7. Control of CHX Impurities

CQA Parameter	Specification for control
pH	<ul style="list-style-type: none"> Drug product at release is controlled to pH 5.5–6.5.
4CA	<ul style="list-style-type: none"> Drug product at release contains <150 ppm.

Aid for CTD Module 3.2.P.3.4 – Analytical and Specifications for CHX Gel

Table 8. Tests and Specifications

CQA Parameter	In-process Specification	Release Specification	Shelf-life specification (acceptance criteria)
Description		A colorless to yellow, translucent gel essentially free from visible particles	A colorless to yellow, translucent gel essentially free from visible particles
Identification of CHX digluconate by: HPLC Ultraviolet		<p>a. Retention time of CHX peak in the sample solution should be within $\pm 3\%$ of retention time of CHX in reference solution preparation.</p> <p>b. Ultraviolet maxima of CHX should be within ± 3 nm relative to reference standard.</p>	
pH		5.5–6.5	5.5–7.0
Apparent viscosity by viscometer (cPs)		3000–8000 cPs	3000–8000 cPs
Minimum fill	<p>Target: 3.3 g</p> <p>Range: 3.0–3.6 g</p>	As per USP <755> Minimum Fill	
CHX digluconate		6.75–7.45 % w/w	6.39–7.45% w/w

content by HPLC (%w/w)		(95.0%–105.0% of label claim)	(90.0%–105.0% of label claim)
Seal Integrity	Confirm sachet integral (pass/fail)		
Sachet appearance	Confirm satisfactory sachet appearance		

Aid for CTD Module 3.2.P.7 – Manufacturing and Packaging

Table 9. Packaging and Composition

Unit Operation	Comments
Sachet Laminate Structure to Conform with European Directive No. 10/2011 for Food Contact and Ph. Eur. 3.3.3 Polyolefins	Outside Layer
	Polyethylene terephthalate (12 µm)
	Low-density polyethylene white pigmented extrusion laminate
	Low-density polyethylene aluminum foil (9 µm)
	Low-density polyethylene-extrusion laminate
	Low-density polyethylene film (40 µm)
	Product Contact Layer
	Example GSK sachet laminate structure
Extractables/Leachables	<ul style="list-style-type: none"> Undergo a risk assessment: foil laminate material, filling and sealing processing of foil laminate, and gel bulk manufacture.

References

Web site

usp-pqm.org. Promoting Quality Medicines Technology Transfer Report-Chlorhexidine digluconate (7.1%) gel [May 2019] Available from: <https://www.usp-pqm.org/sites/default/files/pqms/article/gsk-chx-gel-tech-transfer-report-6-20-2019.pdf>