Insights Into Cannabis Products Quality

Mahmoud A. ElSohly, PhD, BCFM, BCFE

National Center for Natural Products Research and Department of Pharmaceutics and Drug Delivery, School of Pharmacy, University of Mississippi, University, MS 38677





NIDA Marijuana Project Activities

- Grow, harvest, and process cannabis plant material to produce standardized marijuana of different potencies for research
- Isolate and characterize different cannabis components for pharmacological studies
- Prepare bulk quantities of extracts and specific purified cannabinoids (for example, THC, CBD, CBN, CBC, and CBG)
- GMP preparation of materials for clinical trials
- Analysis of confiscated Cannabis materials

















Cannabis Useable Biomass in Barrels

A : THC Derivative Crystals and B : THC from Cannabis







Cannabis Cigarettes

Cannabis Extract



Quality Control: Cannabis Cigarettes

Preparation of Can	nabis Cigarettes (#11	04)	Lot:	
MPR 3.1.1	Revision 1	TRUE COPY by	г.	Date:
	P Cannabi	reparatio is Cigar	on of ettes #1	104
	Pr	oduction R	ecords	
	LOT	:		
	-			
	Master Production Rec	ord Approved By		Date
Project Director:				
1	Master Production Reco	ord Reviewed By		Date
Project Co-Director	C			
Administrative Coo	rdinator:			
Quality Control Uni	t:			

High THC Cannabis Extract





Quality Control: Cannabis Extract

Preparati	on of Cannabis Extra	ct #1101	Lot:		
- MPR 2.1	Revision 10	TRUE COPY	by:	Date:	
	Р	reparati	ion of		
	Cannak	ois Ext	ract # 1	101	
	1 kg	g to 20 k	g Batch		
			-	,	
	Pro	duction R	ecord 2.1		
		Revision	n 10		
	Lot:	1101-			
Ма	ster Production Rec	cord Approved E	Ву	Dat	e
roject Director:					
Ma	ster Production Red	cord Reviewed E	Зу	Dat	e

CoA: Cannabis Extract

Road ller Complex 48	MISSISSIPI			
ty of Mississippi IS 38677-1848 562-915-5928	School of Pharmacy			
CERTIFIC Cannabis Extract : 1102 1102-1812-01	CATE OF ANALYSIS			
TEST Potency (by GC)	SPECIFICATION	RESULT		
Cannabidiol (CBD)	≥25% w/w	54.5%		
Other Cannabinoids (by GC)				
Tetrahydrocannabivarin (THCV)	≤0.3% w/w	<0.025%*		
Δ9-tetrahydrocannabinol (THC)	≤5% w/w	1.60%		
Cannabinol (CBN)	≤5% w/w	0.30%		
Cannabichromene (CBC)	≤5% w/w	1.73%		
Δ ⁸ -tetrahydrocannabinol (Δ ⁸ -THC)		0.68%		
Cannabigerol (CBG)	≤5% w/w	1.08%		
CPD.TUC action	> 20:1	24.1		
10.0258/1 = Constitution Limit	2 20.1	34.1		
Degradation products (by GC)	_			
Campabinal (CBN)	< 20% what	ΝΔ		
NA = Not applicable	3270 W/W			
Ternenes (by GC)				
Total Terpenes	≤1% w/w	1.24%		
		1. S. S. S.		
Residual Solvent (by GC)	and the second s	and the second second second		
Hexane	≥1% w/w	0.08%		
Less en Duden (hu ID rediction)				
Total Maistura & Valatila Impurities	< 60/ w/w	0.60%		
Total Moisture & Volatile Impurfues	≥0% W/W	0.03%		
Heavy Metals (by ICP-MS)		N TO THE OWNER OF THE OWNER		
Lead (Pb)	≤0.5 ppm	0.100 ppm		
Mercury (Hg)	≤1.5 ppm	0.013 ppm		
Cadmium (Cd)	≤0.5 ppm	0.014 ppm		
Arsenic (As)	≤1.5 ppm	0.014 ppm		
Aflatovina	and had been	Participant and a second		
	<5 nnh	Not Detected		
AFB1, AFB2, AGF1, AFG2	≤20 ppb	Not Detected		
		and a second second		
Microbial	L. L. B. L. B.			
Escherichia coli	Absent	Absent		
Salmonella	Absent	Absent		
Total Aerobic Microbial Count (TAMC) <2021>	≤10 ³ cfu/g	< 10 cfu/g		
Total Yeast & Mold Count (TYMC) <2021>	≤10 ² cfu/g	< 10 cfu /g		
	1	28		
Store in freezer.				
Expires: 12/2023				

Quality Control: GMP CBD

MPR 2.3	Revision 2	TRUE COPY hur	Date	y.
	1 to Aldron is	the out of.	Jan	**
	P	reparation o	f	
	Canı	nabidiol # [·]	1003	
	0.1 gram	to 500 gram	is Batch	
	Pro	duction Record Revision 2	2.3	
	PRODUCT			
	NAME: CA	ANNABIDIOL (CBD)		
	~			
Ma	ster Production Re	cord Approved By		Date
ect Director:				
Ma	ster Production Re	cord Reviewed By		Date

Quality Control: GMP THC

Preparation of	f ƻ-tetrahydrocanna	abinol #1001	Lot: 1001 -	-
- MPR 3.1	Revision 1	TRUE COPY by		Date:
	F	Preparatio	on of	
Δ ⁹ -	tetrahy	drocanr	abinol	#1001
	1 gram t	to 1000 g	rams Bat	ch
	-			
	Pro	duction Re Revision	cord 3.1 1	
	Lot:	1001 -		
	-			
	ster Production Re	cord Approved By	1	Date
Mas				
Mas Project Director:				

Development and Validation of a GC-FID Method for the Quantitation of 20 Different Acidic and Neutral Cannabinoids



GC/FID Conditions

- DB-1MS column (15 m x 0.25 mm, and 0.25 µm film thickness).
- Helium (carrier gas) at a flow rate of **0.8 mL/**min.
- The inlet temperature was set at 275°C with a split ratio of 20:1.
- The temperature program started at **190°C to 300°C**
- The total run time was 17.5 min. GC/FID Conditions



Sample Preparation

- Triplicates of the ground samples (100 mg each)
- Extraction with 10 mL of a MeCN :Me OH mixture (8:2) by sonication for 30 min.
- The mixture **was centrifu**ged for **5 min**. at 1,252 xg and transferred into pre-labeled extraction tubes. Aliquots of 10 μ L, 50 μ L, and 100 μ L --- to GC vials for silvation.
- A solution of 4-androstene-3,17-dione at 1 mg/mL was prepared in methanol—chloroform (9:1, v/v) is used as Internal standard .





Silylation of the Samples after extraction

- To each vial, 50 μL of 1 mg/mL I.S. solution and 10 μL of 2% DMAP were added and the solvents evaporated to dryness using a gentle flow of nitrogen gas at 50° C.
- The residue was then silvlated by adding 100 μL of BSTFA, vortexed, and the capped vials were kept in a 70
 ° C oven for 30 min.
- The vials were then brought to room temperature and the contents were transferred to 100 μL GC-vial inserts and analyzed by the GC-FID.







Chemical structures of the 20 silylated acidic and neutral cannabinoids.



Standard 20 silylated cannabinoid mixture (25 µg/mL) and IS (50 µg/mL)

pA SWL-ACCO SWL-ACCO SUL-ACCO	4.833 - THCV-TMS	-5.401 - CBT-TMS 	6:064 - CBC-TMS 6:243 - D8-THC-TMS 6:408 ∆9-THC-TMS	6.718 - CBDV-TMS 7.024 - CBN-TMS 7.024 - CBN-TMS 7.406 - CBC-TMS	7.877 - THCVA-TMS 8.158 - CBDA-TMS 8.383 - CBLA	S 133 133 133 133 133 133 133 133 133 13	9.623 - THCAA-TMS 9.837 - CBCA-TMS 10.222 - CBGA-TMS 10.222 - CBGA-TMS 10.428 - CBEA-TMS
4	5		6	7	8	9	10 min



Cannabinoids Concentration (High CBD Chemovar)



Cannabinoids Concentration (High THC Chemovar)





Cannabinoids Concentration (Intermediate Chemovar)

	FID1 A, FID1A, front detector (F//GC D DR18 HIGH THC_THC_IM_VARIETIES 2022-08-18 15-54-32/046F2601.D)		
Norm]	S S		
225		SMT	
200		САА.	
175 -		Ē	
150		.644 MS	SMS M
125 -		6 II-W	GA-1
100		CBC	B B
75	8.74; (2.86 C. C. C. C. 2. 74 · C. 2. 17	38 -	221-
50		9.8	01 01
25			<u> </u>
4		10	min



Development and Validation of a GC-FID Method for the Quantitation of Δ⁸-tetrahydrocannabinol and other Cannabinoids present in commercial **Δ⁸-tetrahydrocannabinol**



Chemical structures of the identified compounds



OH



OH

Olivetol (1)

HO

iso-THCBF (5)

Ô







 9^{α} -hydroxy-HHC (**13**)

H

9β-hydroxy-HHC (14)





CBT (2)

-0 H 11 Δ^8 -cis-iso-THC(3)

GC-FID representative chromatogram of the 14 cannabinoids at 100 µg/mL and the IS at 10 µg/mL





Standard 20 silylated cannabinoid mixture (25 µg/mL) and IS (50 µg/mL)

- Twenty-one Δ^{8} -THC vape samples were analyzed
- From each Δ^8 -THC vape product, 50 mg was weighed
- Dissolved in 5 mL methanol, vortexed for 10 seconds, sonicated for 5 minutes, then transferred to a 10 mL volumetric flask.
- The volume was adjusted to mark with methanol to get a final concentration of 5 mg/mL.
- To each vial was added 10 μL of the IS (1mg/mL) and the volume adjusted to 100 μL .
- Anthracene at 1 mg/mL was used as IS



GC-FID chromatogram of a representative sample prepared at 5 mg/mL (sample # EA 324)





Acknowledgements

• My Co workers-

Mohamed Radwan, Amira Wanas, Ikhlas Khan, Soumyajit Majumdar Waseem Gul, Don Stanford, Chandrani Gon, Suman Chandra, Hemant Lata.

- The support of School of Pharmacy administration, Business office and the NCNPR Director's office is greatly appreciated.
- Supported in part by National Institute on Drug Abuse (NIDA) contract # NO1DA-10-7773.







